

Re-useable GRAY PAPR Hood

(Multiple staff, multiple patients – Not disposable)

Assembly




Step 1



Step 2



Step 3



Step 4

Step 1

Adjust helmet to your head

Step 2

Place helmet inside re-useable PAPR unit

Step 3

Pull the air inlet through the back of the hood. Make sure it is BELOW the blue airflow adjustment (in the groove)

Step 4

Attach blue pinch clip onto the breathing tube

*Perform hand hygiene before & after donning and doffing your PAPR

Cleaning

- A clean cloth or sponge dampened with a mild solution of water and liquid soap may be used to wipe down the hoods, headcovers, plastic head suspensions, and other plastic parts.
- Hang hoods and headcover by the loop sewn into the top (See image for Step 3) or lay flat.
- Never fold or crush the visor.
- Air dry inside and out thoroughly before storage or reuse.



Hover over the QR code for a YouTube demonstration of the assembly. Start at 2:05 for this type of hood.

Single Person WHITE PAPR Hood

(ONE Staff person, Multiple patients, Multiple shifts)



Donning

- Fasten belt of PAPR unit around waist, draping hose over your shoulder
- Connect hose to back side of PAPR and turn on
- Pull PAPR on over your head. Follow PAPR Guidelines found in the Respiratory Protection Program course in NM Learning & Performance

*Perform hand hygiene before & after donning and doffing your PAPR

Doffing & Cleaning

- Carefully remove PAPR hood in Ante Room, or outside of patient room
- Turn off PAPR unit
- Disconnect hose from back of hood
- After each use, the PAPR hood should be wiped down with purple top Sani Wipes & allow to dry
- Store in a brown paper bag (with your name on it) between uses in a designated areas on your unit.
- Hood should be discarded & replaced if any damage is detected

Process for Obtaining and Returning PAPR (3 South Pyxis)

Purpose: to assure we have PAPRs available for those who require them

S: PAPRs cannot be easily located when removed under temporary patient/floorstock

B: Process to get a PAPR is listed below

A: The patient is not being selected when retrieving a PAPR

R: Always **choose the patient** that requires the use of a PAPR



PAPR TAKE/RETURN

Get a PAPR

1. Log into Pyxis
2. Select correct patient
3. Choose **Select Items**
4. Take 1 blue PAPR Bag
 - a. press **TAKE** button once
5. Take 2 hoods
 - a. press **TAKE** button twice
6. Log out of Pyxis
7. Choose a battery that matches the number on the PAPR bag

BATTERIES/CHARGERS ARE TO THE LEFT, NEXT TO THE FRIDGE

Return a PAPR

1. Log into Pyxis
2. Select correct patient
3. Choose **Select Items**
4. Return 1 blue PAPR Bag into the correct space, matching the number on the bag to the Pyxis space
 - a. press **RETURN** once
5. Log out of Pyxis
6. Return the battery that matches the number on the PAPR bag

BATTERIES/CHARGERS ARE TO THE LEFT, NEXT TO THE FRIDGE