





North Region N95 Selection

Note: This information is subject to change based on inventory levels

Background: N95 mask selection should ensure appropriate seal for each individual while conserving supply. To ensure availability of appropriate masks for all, please follow the below guide to select an appropriate N95 Mask

DESCRIPTION	IMAGE	INTENDED USER
Manufacturer: 3M Model: 8210 Size: Regular		All staff/providers who are able to achieve a proper seal ¹ . Note: This mask must be worn with a face shield for splash protection
Manufacturer: Moldex Model: 1511, 1512, 1513 Size: Small, Medium, Large		All staff/providers who cannot achieve a proper seal ¹ with the 3M 8210.
Manufacturer: 3M Model: 1870 Size: One size		All staff/providers who cannot achieve a proper seal ¹ with the 3M 8210 or Moldex 1511, 1512, or 1513 or Kimberly Clark/Halyard respirator or the 3M 1860 respirator. <i>* These respirators are currently on backorder with manufacturer with a very limited supply on hand. Manufacturer has no current release date available and will be phased out</i>
Powered Air-Purifying Respirator (PAPR)		All providers/staff who cannot achieve a proper seal with either of the above masks. Note: Reach out to your manager for directions on how to obtain a PAPR

1. Refer to [Seal Check Instructions](#) on [Physician Forum PPE Resources page](#) to ensure proper fit
2. **Note:** If individuals have respirator fit concerns or cannot properly achieve a respirator fit, contact Corporate Health and make an appointment for a respirator fit test.